

# Chromatin Immunoprecipitation (ChIP)

Aidan Quinn

December 22, 2018

## Contents

<b>1</b>	<b>Introduction</b>	<b>1</b>
<b>2</b>	<b>Required Reagents and Equipment</b>	<b>2</b>
2.1	Reagents . . . . .	2
2.2	Equipment . . . . .	2
<b>3</b>	<b>Protocol</b>	<b>3</b>
3.1	Crosslinking . . . . .	3
3.2	Chromatin Sonication . . . . .	3
3.3	Immunoprecipitation . . . . .	4
3.4	Pull-down . . . . .	4
3.5	Reverse crosslinking . . . . .	4
3.6	ChIP DNA Clean-up . . . . .	5
<b>4</b>	<b>Recipes</b>	<b>6</b>
4.1	Covaris Lysis Buffer . . . . .	6
4.2	Covaris Wash Buffer . . . . .	6
4.3	Agilent Sonication Buffer . . . . .	6
4.4	Blocking Buffer . . . . .	6
4.5	Wash Buffer A . . . . .	7
4.6	Wash Buffer B . . . . .	7
4.7	Elution Buffer . . . . .	7

## 1 Introduction

This protocol aims to describe the standard chromatin immunoprecipitation protocol using **Agilent sonication buffer** with sonication in the Covaris sonicator and using **magnetic bead** pull-down. This protocol is standardized for histone and histone modification ChIP, additional considerations may be necessary for application of this protocol to transcription factor or TF co-factor ChIP (*i.e.* consider micrococcal nuclease digestion in place of sonication).

*A note on condition optimization:*

The two most critical parameters to optimize are the sonication and the fixation conditions, with the former predominating. These conditions should be fixed for a given sample type or cell line but may vary quite significantly across cell/tissue types and should be set up carefully.

## 2 Required Reagents and Equipment

### 2.1 Reagents

- Lysis buffer (Covaris)
- Wash buffer (Covaris)
- Wash buffer A
- Wash buffer B
- Sonication buffer (Agilent)
- Blocking buffer
- ChIP IP buffer
- Elution buffer
- Methanol-free formaldehyde
- Glycine (1.25 M)
- PBS
- Dyna Beads (coated with protein A, G or A/G)

### 2.2 Equipment

- Magnetic bead separation rack (pre-cooled at 4 °C)
- Rotating rack stored in 4 °C closet
- Covaris sonication system
- Low-retention pipette tips
- Low-retention 1.5 ml conical tubes

## 3 Protocol

### DAY 1

[*n.b.*] Before beginning, check that all buffers (Covaris lysis, wash, and Agilent sonication) contain fresh protease inhibitor cocktail.

#### 3.1 Crosslinking

1. Harvest cells in log phase growth. Collect a maximum of  $50 \times 10^6$  cells in 30 ml.
2. Add 2 ml of 16% methanol-free formaldehyde (final 1% formaldehyde). Invert several times and incubate 10' (RT) on a shaking platform.
3. Stop the crosslinking reaction by adding glycine at a final concentration of 0.125 M (3.5 ml/tube from 1.25 M glycine stock). Incubate for 5' (RT) on a shaking platform.
4. Centrifuge 5' @ 500 g (4°C), aspirate supernatant and place cells on ice.
5. Wash cells twice in ice-cold PBS with 5'/500 g (4°C) centrifugation.
6. **\*\*Optionally freeze fixed cells\*\*** **Warning:** it is preferable to continue through the sonication and instead freeze sonicated chromatin.
  - (a) Completely aspirate supernatant from fixed cell pellet and snap-freeze in  $LN_2$  or dry-ice ethanol bath.
  - (b) Store fixed cell pellet at -80°C.

#### 3.2 Chromatin Sonication

1. Resuspend fixed cell pellet (max  $100 \times 10^6$  cells) in 10 ml of ice-cold Covaris lysis buffer in a 15 ml conical tube. Incubate rotating 10' (4°C).
2. Centrifuge 5'/4000 rpm (4°C) to pellet nuclei.
3. Wash nuclear pellet (resuspend) in 10 ml ice-cold Covaris wash buffer, incubate rotating for 10' (4°C), then centrifuge 5'/4000 rpm (4°C) to pellet nuclei.
4. Resuspend the nuclear pellet in the appropriate amount of ice-cold Agilent sonication buffer to achieve 10 – 20  $\times 10^6$  cells/ml.
5. Sonication should be performed uniformly (constant volume, cell concentration, cell number, *etc.*) in up to 1 ml total volume per sonication tube using the Covaris sonicator. Sonication time should be carefully optimized to achieve a fragment length distribution between 200-500 bp.
6. **\*\*Optionally freeze chromatin\*\***
  - (a) Split sonicated chromatin into eppendorf tubes and centrifuge 15'/13,000 rpm (4°C)
  - (b) remove supernatant and place in a new eppendorf tube for snap-freezing (in  $LN_2$  or dry ice ethanol bath) and store at -80°C.

### 3.3 Immunoprecipitation

1. Pool sonicated chromatin samples (by biological condition) in a 15 ml conical tube. Add TritonX-100 to 1% final concentration. Briefly vortex and divide equally among eppendorf tubes according to the total number of IP reactions.
2. Set aside 50  $\mu$ l of sonicated chromatin as input and store at 4 °C until IP is complete.
3. Add the appropriate antibody (Ig control or specific Ab). *n.b.* Antibody amount needs to be empirically determined, normally in range 1-6  $\mu$ g per IP.
4. Incubate overnight at 4 °C with rotation.

#### DAY 2

### 3.4 Pull-down

1. Wash 100  $\mu$ l of protein A (or G) magnetic beads three times in 1 ml of ice-cold blocking buffer by inverting several times in a 1.5 ml tube. After each wash, collect beads using a pre-cooled magnetic rack.  
[*n.b.*] If downstream analysis is qPCR, 50  $\mu$ l of beads is sufficient, however for ChIP-seq 100  $\mu$ l is required.
2. Transfer the washed beads to the tube containing chromatin+Ab.
3. Incubate on a rotating platform at 4 °C for 4-5 hours.
4. Remove supernatant from beads using the pre-cooled magnetic stand. Add 1 ml of ice-cold ChIP IP buffer (sonication buffer + 1% TritonX-100) and incubate with rotation 5' (4 °C). Repeat once.
5. Remove supernatant from beads using the pre-cooled magnetic stand. Add 1 ml of ice-cold Wash buffer A and incubate with rotation 5' (4 °C). Repeat once.
6. Remove supernatant from beads using the pre-cooled magnetic stand. Add 1 ml of ice-cold Wash buffer B and incubate with rotation 5' (4 °C). Repeat once.
7. Remove supernatant from beads using the pre-cooled magnetic stand. Add 1 ml of ice-cold 1X TE buffer (pH 8.0) with 50 mM NaCl. Incubate 10' (4 °C) with rotation.
8. Remove supernatant from beads using the pre-cooled magnetic stand. Add 210  $\mu$ l of Elution buffer to the beads. Mix by *delicate* vortexing.
9. Incubate beads in Elution buffer 16' (65 °C) resuspending by *delicate* vortexing every 2'.
10. Centrifuge 2'/10,000 rpm (RT) and recover the supernatant in new, low-binding tubes.  
[!] From this point forward use only **low binding** tips and tubes.
11. Add 150  $\mu$ l of Elution buffer to the input chromatin (from previous day).

### 3.5 Reverse crosslinking

1. Reverse crosslinking in both the IPs and input by incubating in 65 °C water bath over-night (12-14 hours).

**DAY 3****3.6 ChIP DNA Clean-up**

1. Add 200  $\mu$ l of 1X TE buffer (pH 8.0, 50 mM NaCl) to each tube to dilute the SDS.
2. Add 40  $\mu$ l of RNase A (molecular grade) and incubate at at least 2 hours (37°C).
3. Add 80  $\mu$ g of proteinase K (4  $\mu$ l of our stock) and incubate 30' (55°C).
4. Perform DNA clean-up using Qiagen PCR Purification kit.
5. Recover ChIP DNA in low-binding tubes. Recovered ChIP-DNA and Input DNA can be stored at -80°C or immediately used for qPCR or HTS library prep.

## 4 Recipes

\* should be added immediately prior to use.

‡ use only freshly prepared stock.

### 4.1 Covaris Lysis Buffer

	Stock	Final	Amount
Hepes-KOH pH 7.5		50 mM	
NaCl		140 mM	
EDTA pH 8.0		1 mM	
Glycerol		10%	
Triton X-100		0.25%	
NP40*		0.5%	
Protease Inhibitors*		1X	

### 4.2 Covaris Wash Buffer

For 50 mL:

	Stock	Final	Amount
Tris-HCl pH 8.0		10 mM	
NaCl		200 mM	
EDTA pH 8.0		1 mM	
EGTA pH 8.0		0.5 mM	
Protease Inhibitors*		1X	

### 4.3 Agilent Sonication Buffer

For 50 mL:

	Stock	Final	Amount
Tris-HCl pH 8.0		10 mM	
NaCl		100 mM	
EDTA pH 8.0		1 mM	
EGTA pH 8.0		0.5 mM	
N-lauroylsarcosine Sodium Salt		0.5%	
Na-deoxycholate‡*		0.1%	
Protease Inhibitors*		1X	

### 4.4 Blocking Buffer

Filter and store at 4 °C. For 50 mL:

	Stock	Final	Amount
PBS		1X	
BSA		0.5%	

#### 4.5 Wash Buffer A

For 50 mL:

	<b>Stock</b>	<b>Final</b>	<b>Amount</b>
<b>Hepes</b> pH 7.8		50 mM	
<b>NaCl</b>		500 mM	
<b>EDTA</b>		1 mM	
<b>Triton X-100</b>		1%	
<b>Na-deoxycholate</b> ‡*		0.1%	
<b>SDS</b> *		0.1%	
<b>Protease Inhibitors</b> *		1X	

#### 4.6 Wash Buffer B

For 50 mL:

	<b>Stock</b>	<b>Final</b>	<b>Amount</b>
<b>Tris-HCl</b> pH 8.0		20 mM	
<b>LiCl</b>		250 mM	
<b>EDTA</b>		1 mM	
<b>NP40</b>		0.5%	
<b>Na-deoxycholate</b> ‡*		0.1%	
<b>Protease Inhibitors</b> *		1X	

#### 4.7 Elution Buffer

For 50 mL:

	<b>Stock</b>	<b>Final</b>	<b>Amount</b>
<b>Tris-HCl</b> pH 8.0		50 mM	
<b>EDTA</b>		10 mM	
<b>SDS</b> *		1%	